Auto-tRNAseq: Driving throughput in tRNA library generation with firefly[®]

Taylor C. Fennelly¹, Russel M. Vincent^{2,3}, Huw Rees⁴, Russell Buckley-Taylor⁴, Robert Steen¹, and George Church^{2,3} ¹Biopolymers Facility, Harvard Medical School, Boston, USA, ²Department of Genetics, Harvard Medical School, Boston, MA 02115, USA, ³Wyss Institute for Biologically Inspired Engineering, Harvard University, Boston, MA, USA, ⁴SPT Labtech, Melbourn, Cambridgeshire, UK

Introduction

Transfer RNAs (tRNAs) are the critical adaptor molecules that translate genetic information into functional protein sequences. While the genome defines a static set of tRNA gene sequences, the functional tRNA pool depends on the expression and aminoacylation levels of different tRNA species, as well as many modifications that fine-tune tRNA activity.^{1,2} Dynamically regulated tRNA pools dictate the fidelity of the translation process and proteome integrity of a cell. Considering the central role of tRNA in translation regulation, investigating the cellular composition of tRNA pools provides a deeper understanding of the translational state of the cell. High-resolution sequencing-based methods for tRNA sequencing have provided transformative insights into the physiological impact of tRNA expression, modification, and aminoacylation levels on cellular stress and/or disease conditions.³⁻⁶ More recently, large-scale genome engineering efforts have also employed tRNA sequencing to recapitulate the translational state of engineered cells.⁷

While tRNA sequencing offers comprehensive insights into the cellular state, it has not been readily employed in research across the globe. Current library generation workflows for tRNA sequencing are time-consuming, arduous, and rely on multiple gel purification steps that incorporate experimental variability and limit experimental throughput. In this poster, for the first time we propose a high-throughput automated library generation workflow for tRNA sequencing using SPT Labtech's firefly liquid handling platform.

Automating genomics liquid handling with firefly

firefly has been designed specifically to streamline NGS library preparation by bringing together multiple liquid handling capabilities.



16 total deck positions across two moving decks

Additional functionalities (not shown) Below-deck heating and shaking modules Gripper for hands-free plate transfer

sptlabtech.com



Figure 1. Manual workflow for generation of tRNA sequencing data. Chemical treatment steps and bead-based cleanups result in a highly specialized and labor-intensive workflow.

Automated tRNA library generation



Reagent reservoirs

Figure 2. Automated workflow for generation of tRNA sequencing data. Chemical treatment steps and bead-based cleanups are performed on an automated liquid-handler, improving consistency and allowing for higher-throughput experiments

Bead-based RNA purification

firefly was used to perform all dispense, aspiration and mixing steps in a 96-well PCR plate (Biorad Hard Shell Plate). The non-contact dispense head was used to dispense RNA purification beads (SPRIselect, Beckman Coulter), 20-45% isopropanol or ethanol in water and elution buffer (RNase-free water). The pipetting head was used to perform the initial RNA transfer, all mixing steps, and all transfers of supernatant to waste. A 96-well plate magnet (MagnumFLX, Alpaqua) was used for all bead separation steps.

AutotRNAseq_worflow1_part1: (a) RNA size selection: Using the SPRI beads with 7.5% PEG + 20% isopropanol to remove large RNA molecules > Increase isopropanol concentration (~45%) to bind small RNA molecules > (b) Chemical treatment: Add chemical reaction buffer I > Incubate for 15 mins > Add reaction buffer II.

> Dispense isopropanol, RNase-free water, and SPRI beads into large tRNA

Transfer supernatant with size-selected tRNA (small tRNA)



Figure 3. Execution of tRNA size-selection on firefly

Results Small RNA enrichment on firefly





Figure 4. TapeStation results showed successful and consistent size selection of tRNA between all wells in a column. Larger RNA species were removed from the mixture to leave behind only the smaller RNA species of interest.

Conclusions

- firefly can successfully be used to automate bead-based RNA purification, a key component of Next Generation Sequencing (NGS) workflows, with a uniform fragment size distribution demonstrated for samples across a 96-well plate.
- firefly completed size selection of RNA for the purposes of tRNA sequencing and modification detection.

The versatility of firefly with its two distinct liquid handling heads is well-suited to develop customized bead-based RNA size selection methods to circumvent gel purification techniques. In this work, we have prototyped the crucial short RNA (<200 nt) enrichment step with consistent final yield from the same starting material of total RNA input, with minimal well-to-well variation.

References

- 1. Bednářová, A. et al. Lost in Translation: Defects in Transfer RNA Modifications and Neurological Disorders. Front. Mol. Neurosci. 10, (2017).
- 2. Suzuki, T. The expanding world of tRNA modifications and their disease relevance. Nat. Rev. Mol. Cell Biol. 22, 375–392 (2021).
- 3. Zheng, G. et al. Efficient and quantitative high-throughput tRNA sequencing. Nat. Methods 12, 835–837 (2015).
- 4. Evans, M. E., Clark, W. C., Zheng, G. & Pan, T. Determination of tRNA aminoacylation levels by high-throughput sequencing. Nucleic Acids Res. 45, e133 (2017).
- 5. Behrens, A., Rodschinka, G. & Nedialkova, D. D. High-resolution quantitative profiling of tRNA abundance and modification status in eukaryotes by mimtRNAseq. Mol. Cell 81, 1802-1815.e7 (2021)
- 6. Watkins, C. P., Zhang, W., Wylder, A. C., Katanski, C. D. & Pan, T. A multiplex platform for small RNA sequencing elucidates multifaceted tRNA stress response and translational regulation. Nat. Commun. 13, 2491 (2022).
- 7. Schindler, D. et al. Design, construction, and functional characterization of a tRNA neochromosome in yeast. Cell 186, 5237-5253.e22 (2023).

sptlabtech

Fraction of large RNA (>200 nt)

Fraction of small RNA (<200 nt)