

# miniaturised magnetic bead clean-ups on mosquito<sup>®</sup> HV genomics

- new magnetic bead separation block and protocol optimized for bead clean-ups with mosquito<sup>®</sup> HV genomics
- miniaturize NGS library preparation workflows with significant cost savings

An essential accessory for every NGS lab, the magnetic bead separation block has been optimized for miniaturized SPRI bead clean-ups with mosquito HV genomics. The block and the corresponding mosquito protocol provide a semi-automated process that combines significant cost savings with high recovery rates.

The new bead clean-up process will be helpful in single cell and bulk library preparation workflows, including the latest on-bead tagmentation methods.

The improved magnetic bead clean-up process is compatible with all mosquito Genomics models. The block fits also **dragonfly<sup>®</sup> discovery**. It delivers high DNA recovery rates and reliable clean-up or size selection in 384-well PCR plates (Table 2 and Fig.2). Bead cleanups have been optimized to volumes as low as 5  $\mu$ L.

## technical specification


### Tested plates:

Eppendorf twin.tec 384-well PCR plate  
Bio-Rad Hard-Shell 384-well PCR plate

### Table 1:

Recommended volumes for miniaturized bead clean-up process performed with SPT Labtech magnetic bead separation block.

Sample volume	Bead volume	EtOH volume	Elution volume
2-10 $\mu$ L	2-10 $\mu$ L	4.5-10 $\mu$ L	2-20 $\mu$ L




**Fig. 1:** Separation of magnetic beads in a 384-well PCR plate. 5  $\mu$ L of SPRIselect beads were added to 5  $\mu$ L of sample and mixed thoroughly before placing on the magnetic bead separation block (SPT Labtech).

**A:** Immediately after resuspension.

**B:** 30 seconds after the separation process.

### validated with:

- **AMPure XP beads**  
(Beckman Coulter)
- **SPRIselect beads**  
(Beckman Coulter)
- **Illumina DNA Prep**  
(previously known as Nextera DNA Flex)



Purification 0.5x	Full volume manual	Miniaturized manual	Miniaturized mosquito HV
Recovery (%)	67 +/- 1.0	62 +/- 0.7	69 +/- 0.8
Clean-up rate (%)	100	100	100

Purification 0.8x	Full volume manual	Miniaturized manual	Miniaturized mosquito HV
Recovery (%)	72 +/- 0.3	70 +/- 0.8	76 +/- 0.3
Clean-up rate (%)	100	100	100


Size selection	Full volume manual	Miniaturized manual	Miniaturized mosquito HV
Recovery (%)	50 +/- 0.3	45 +/- 0.5	61 +/- 0.9
Clean-up rate (%)	100	100	100

**Table 2:** Control DNA (Mid Range DNA Ladder, Jena Bioscience # M-203, 2.6 ng/μl) was treated with SPRIselect beads according to the manufacturer's Purification or Size Selection protocols (Beckman Coulter). Bead ratios (-fold volume of bead suspension vs. initial sample volume) used: Purification 0.5x or 0.8x, respectively; Size Selection: first step 0.95x; second step 0.85x. Samples were analysed on a Fragment Analyzer with HS NGS Fragment Kit (Advanced Analytical Technologies).

**Full-volume manual:** 50 μl of control DNA in 500 μl PCR-grade microtubes.

**Miniaturized manual:** 5 μl of control DNA in Eppendorf twin.tec 384-well PCR plate with magnetic bead separation block (SPT Labtech).

**Miniaturized, automated:** as in miniaturised manual process but all pipetting steps carried out on mosquito HV genomics.



**Fig. 2:** Control DNA (Mid Range DNA Ladder, Jena Bioscience # M-203, 2.6 ng/μl) was treated according to the SPRIselect purification protocol with a 0.8x volume of bead suspension. Samples were analysed on a Fragment Analyzer using the HS NGS Fragment Kit. Representative results are shown.

**A: Input control** (before clean-up)

**B: Full-volume manual** purification with 50 μl of control DNA in 500 μl PCR-grade microtubes.

**C: Miniaturized manual** purification with 5 μl of control DNA in Eppendorf twin.tec 384-well plate with magnetic bead separation block.

**D: Miniaturized, automated** purification: as in C but all pipetting steps carried out on mosquito HV genomics.

To order, email [consumables@sptlabtech.com](mailto:consumables@sptlabtech.com). Part Number: 3268-02008

**sptlabtech.com**