

comPACT®

automated cold sample storage
and retrieval has never been
more accessible



Scan me to
find out more
about comPACT



Function	comPACT
Temperature	-20°C
Dimension (W X D X H)	1.35 x 0.8 x 1.98 m (53 x 32 x 78 inches)
Capacity	40,000 x 1.4 mL vials
Throughput	14 racks in and out per day (total 2,688 vials per day)

Proven technology in a smaller package

Robust and user-friendly storage in a more affordable design, meaning more labs can benefit from automated sample management.

Assured sample integrity

Tube transport harnesses the power of pneumatics for secure operation. The ability to cherry pick samples avoids unnecessary harmful freeze/thaw cycles.

Grows with your business needs

Start small now and expand later with comPACT's flexible and scalable modular structure. Modules can be located anywhere and samples delivered remotely.

Complete operational efficiency

Simple interface connects easily to LIMS and automated sample processing workflows.



comPACT's vending machine-style retrieval delivers chosen samples directly to you without the stress of manual searching